

Monarch[®] Plasmid Miniprep Kit Protocol Card




NEB #T1010

For a detailed protocol or to download the full manual, visit www.neb.com/T1010.

BEFORE YOU BEGIN:

- Add 4 volumes of ethanol ($\geq 95\%$) to one volume of Plasmid Wash Buffer 2.
- All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).
- If precipitate has formed in Lysis Buffer (B2), incubate at 30–37°C, inverting periodically to dissolve.
- Store Plasmid Neutralization Buffer (B3) at 4°C after opening, as it contains RNase A.
- If working with plasmids ≥ 10 kb, preheat the appropriate amount of DNA Elution Buffer to 50°C.

PROTOCOL STEPS:

- 1. Pellet 1–5 ml (not to exceed 15 OD units) bacterial culture by centrifugation for 30 seconds. Discard supernatant.** 1.5 ml of culture is sufficient for most applications. Ensure cultures are not overgrown (12-16 hours is ideal).
- 2. Resuspend pellet in 200 μ l Plasmid Resuspension Buffer (B1)** . Vortex or pipet to ensure cells are completely resuspended. There should be no visible clumps.
- 3. Add 200 μ l Plasmid Lysis Buffer (B2)** , gently invert tube 5–6 times, and incubate at room temperature for 1 minute. Color should change to dark pink, and solution will become transparent and viscous. Do not vortex.
- 4. Add 400 μ l of Plasmid Neutralization Buffer (B3)** , gently invert tube until neutralized, and incubate at room temperature for 2 minutes. Sample is neutralized when color is uniformly yellow and precipitate forms. Do not vortex.
- 5. Centrifuge lysate for 2–5 minutes.** For best results, and especially for culture volumes > 1 ml, we recommend a 5 minute spin to ensure efficient RNA removal by RNase A. Pellet should be compact; spin longer if needed.

continued on back →

6. Carefully transfer supernatant to the spin column and centrifuge for 1 minute. Discard flow-through.
7. Re-insert column in the collection tube and add 200 μ l of Plasmid Wash Buffer 1. Centrifuge for 1 minute. Discarding the flow-through is optional.
8. Add 400 μ l of Plasmid Wash Buffer 2 and centrifuge for 1 minute.
9. Transfer column to a clean 1.5 ml microfuge tube. Use care to ensure that the tip of the column does not come into contact with the flow-through. If there is any doubt, re-spin the column for 1 minute.
10. Add ≥ 30 μ l DNA Elution Buffer to the center of the matrix. Wait for 1 minute, then spin for 1 minute to elute DNA. Nuclease-free water (pH 7-8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA, (≥ 10 kb), heating the elution buffer to 50°C prior to use can improve yield.

Questions?

Our tech support scientists would be happy to help.

Email us at info@neb.com

This product is covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information about commercial rights, please email us at gbd@neb.com. The use of this product may require you to obtain additional third party intellectual property rights for certain applications.

New England Biolabs is an ISO 9001, ISO 14001 and ISO 13485 certified facility.

© Copyright 2019, New England Biolabs, Inc.; all rights reserved.

This card is made with FSC certified 100% post-consumer fiber. Please recycle.

V3.0 - 7.19 #101025

