


# TOXOPLASMOSIS (2)


Indirect Immunofluorescence (IFI)  
Titration of anti-*T.gondii* IgG



*Toxoplasma* = Ag



# Evolution kinetics of anti-*Toxoplasma* antibodies




Seroconversion = négative serology → positive

The appearance of IgG allows to affirm the seroconversion

# Slides used for IFI

- Antigen bound to the slide :  
*Toxoplasma gondii*
- Slide of 5 wells: 1 for 2 groups



# Dilution of serums

- The serums Tx1 to Tx6 are given diluted at 1/4 and have to be diluted at 1/64 in PBS:
- Dilutions are performed in PBS, in a microtitration plate:

–  $\frac{1}{2}$        $\Rightarrow$  **1/8**      50 $\mu$ L+ 50

–  $\frac{1}{2}$        $\Rightarrow$  **1/16**      50      + 50

–  $\frac{1}{2}$        $\Rightarrow$  **1/32**      ...

–  $\frac{1}{2}$        $\Rightarrow$  **1/64**

- Control serums TxT+, TxT-, TRf are ready to use  
→ ask for the deposition of TxT+ on the slide

PBS

50 50 50 50  
μl μl μl μl

[illegible]

## Dilution

1/8

Serum  
(1/4)

50  
μl

PBS

50  
μl

50  
μl

50  
μl

50  
μl

1

2

3

4

5

6

7

8

9

10

1

12

1

A

B

C

D

E

F

## G

H

## Dilution

$$\frac{1}{8} \quad \frac{1}{16}$$

Serum  
(1/4)

50 50  
μl μl

PBS

50      50      50      50  
μl      μl      μl      μl

[illegible]



## Dilution

$$\frac{1}{8} \quad \frac{1}{16} \quad \frac{1}{32}$$

Serum  
(1/4)

50      50      50  
μl      μl      μl

PBS

50      50      50      50  
μl      μl      μl      μl

[illegible]

## Dilution

1/8    1/16   1/32   1/64

Serum  
(1/4)

50      50      50      50  
μl      μl      μl      μl


PBS

50      50      50      50  
μl      μl      μl      μl

[illegible]

# Serum loading

- Put the slide in a humid chamber
- **20 $\mu$ L of diluted serums in the wells as illustrated:**



TxT+ = positive control

TxT- = negative control


TRf = Control reagent

# Incubation

- In a humid chamber
- Incubator **37° C**
- **30 min**

# Washes

- Rinse gently with **PBS** (pissette)



- 1<sup>st</sup> wash in **PBS** : 5 min
- 2<sup>nd</sup> wash in **PBS** : 5 min
- 3<sup>rd</sup> wash in **PBS** : 5 min



# Slides drying

- Very gently !!!
- Between the wells  
as illustrated →

Coton bud

Paper towel



# Secondary antibody loading

- Human anti-IgG conjugated with fluorescein.
- Put the slide in a humid chamber
- Anti-IgG is ready to use: ask for the deposition on the slide

# Incubation

- 37° C for 30 min in a humid chamber

# Washes

- 3 washes 5 min in PBS (as the first time)
- **At the last washing step:** add 3 to 5 drops of Evans Blue in PBS bath

# Drying

- Always gently!!!
- More gently than the first time





# Mounting

- In buffered glycerol
- 5 drops deposited in each well
- Recover with large coverslip 24 x 60 mm
- Remove the excess of mounting medium if necessary
- Possibility to keep the slides mounted several days at +4° C

# Reading

- **Epifluorescence** microscope
- Objective x **40**
- **The interpretation is not possible if the controls are not correct**
- What can we see?



# Analysis of results

- **Membrane fluorescence: positive**
- If fluorescence only in the cytoplasm, or no fluorescence: **negative**
- **Positivity threshold of the technique: 1/64**
- **Conclude** and indicate the **patient follow-up**